

zapper, mora, metatron, rayometer, non-line health system

I'll treat these as terms/devices you want identified and compared, because several are niche and easy to mix up. I'll verify what each one usually refers to, then separate "claimed use" from "scientific/medical evidence."

Thought for 11s >

These terms are mostly from **alternative / energy-medicine / bioresonance** devices, not standard hospital diagnostic systems.

Term	What it usually means	Claimed purpose	Evidence / caution
Zapper	Small electrical device popularized by Hulda Clark	Claims to kill parasites, bacteria, viruses, or support disease treatment	Strong caution. The FTC charged sellers of "Super products for unsupported claims about cancer/A other serious diseases. Fede... +1
MORA	A type of bioresonance therapy device/system, often called MORA therapy	Claims to detect or treat imbalance using body "frequencies"	Some studies exist, but many are retrospective/n controlled or have conflicts. Major cancer center describes bioresonance as unproved for diagnosis/treatment. ScienceDirect +1
Metatron / NLS	"Non-linear diagnostic system" often marketed as a scan/analysis device	Claims to scan organs and detect health problems through frequency/information fields	Be careful. "Metatron" search results are mixed, a health-device claims are vendor/clinic claims, not medical diagnostics. It should not replace real tes blood tests, ultrasound, CT, MRI, ECG, etc.
Rayometer / Rayonex / Rayocomp	Usually connected to Rayonex bioresonance systems, "Bioresonance according to Paul Schmidt"	Claims to harmonize or support self-regulation using frequency spectra	Mostly alternative medicine. Vendor sites describ wellness/self-regulation claims, but this is not the clinically validated diagnosis/treatment. Rayonex B

Term	What it usually means	Claimed purpose	Evidence / caution
Non-linear health system / NLS	Broad name for devices claiming non-linear energetic/body-field analysis	Claims to analyze organs, cells, stress, toxins, parasites, etc.	Treat as non-standard. For real disease diagnosis, licensed medical devices and physicians.

Simple explanation:

These systems are usually in the same family: **bioresonance / frequency / energy-based alternative health devices**. They often claim that the body has electromagnetic or informational frequencies, and the device can detect or correct abnormal patterns.

My recommendation:

Use them only as **wellness/alternative tools**, not as medical diagnostic equipment. Do **not** use them to decide whether you have cancer, infection, organ disease, parasites, diabetes, heart disease, or allergy. For those, use standard medical tests.

A good safety rule: if a device claims it can diagnose many diseases quickly, scan the whole body, or replace lab/imaging tests, be very skeptical.

    ...    Sources

I've heard about zapper, mora, metadron, rayometer, non-line health system in bioresonance field.

Yes, these names are commonly grouped around the **bioresonance / frequency medicine / NLS** market. I'll sort them by category so it's easier to understand what is a device type, what is a brand/system, and what is mainly a marketing term.

Thought for 5s

Yes — in the **bioresonance field**, these names are usually connected, but they are not exactly the same thing.

Name	Category	Meaning in bioresonance field
Zapper	Electrical "frequency" device	Usually a small device that applies weak electrical pulses through hand electrodes or pads. Often linked to Hulda Clark-style claims about parasites, infections, or detox.

Name	Category	Meaning in bioresonance field
MORA	Bioresonance system / therapy method	One of the older European bioresonance systems. "MORA" is commonly said to come from MO rell and RA sche, the names associated with its development. It is often marketed for allergy, pain, chronic conditions, and "body frequency balancing."
Metatron / Metadron	NLS diagnostic system brand/type	Usually refers to Metatron Hunter / Metatron NLS devices. These are marketed as "non-linear diagnostic systems" that scan the body using headphones or sensors and software images.
Rayometer / Rayonex	Bioresonance brand/system	Connected with Rayonex and "Bioresonance according to Paul Schmidt." Often marketed for frequency testing, harmonization, and energetic regulation.
Non-linear health system / NLS	Device category / marketing term	"NLS" usually means Non-Linear System or Non-Linear Screening/Diagnostics . It is used for devices like Metatron that claim to analyze organs, tissues, stress, or disease tendencies.

The common idea behind all of them is this:

The body is claimed to have electromagnetic or informational frequencies, and disease or imbalance is claimed to create abnormal frequency patterns. The device claims to detect, invert, cancel, or harmonize those patterns.

But medically, this is the important part: **bioresonance is not accepted as a reliable diagnostic or treatment method in mainstream medicine**. Memorial Sloan Kettering describes bioresonance therapy as an **unproved** device approach used for diagnosis and treatment claims, and says it is not recommended for preventing or treating cancer.

Memorial Sloan ... +1

For **Zapper** devices specifically, the FTC previously charged sellers connected with Hulda Clark-style "Zapper" products for making unsupported claims, including claims related to serious diseases.

Federal Trade C... +1

For **MORA**, there are some published studies and practitioner reports, but the overall evidence base is weak compared with standard medical diagnostics/treatments. WebMD describes bioresonance as using a device such as the Mora device to measure electromagnetic waves believed to come from the body, while MSK classifies the broader approach as unproved.

WebMD +1

For **Metatron/NLS**, many descriptions come from vendors, not from strong clinical validation. Vendor pages market it as a non-invasive diagnostic scanner based on bioresonance or NLS diagnostics, but that is not the same as being validated like ECG, ultrasound, CT, MRI, blood chemistry, or FDA-cleared diagnostic tools.

Neogenesis Syst... +1

So the clean classification is:

Group	Examples
Simple electrical frequency devices	Zapper
Classical bioresonance therapy systems	MORA, Rayonex/Rayometer
Software-based "body scan" / NLS systems	Metatron, Hunter 4025, Non-linear health system
Broad field name	Bioresonance, frequency medicine, energy medicine, NLS diagnostics

My practical recommendation: treat them as **alternative wellness devices**, not as medical diagnostic equipment. They should not be used to diagnose cancer, infection, parasites, organ disease, allergy, diabetes, heart disease, or replace lab/imaging tests.

    ...  Sources

give me detailed datasheet for each device

I'll make this as a **practical vendor-style datasheet**, but one warning first: these device names are broad families, not exact models. So I'll give a detailed **typical/specification datasheet** for Zapper, MORA, Metatron/NLS, Rayometer/Rayonex, and NLS health systems, and I'll mark where specs depend on the exact model.

Thought for 24s >

Below is a **detailed practical datasheet** for the common device families you mentioned. Because "zapper," "MORA," "Metatron," "Rayometer/Rayonex," and "NLS/non-linear health system" are not one exact model each, I'm giving the **most common representative model/spec family** and marking anything that depends on the supplier.

Important: these are **not accepted replacements for standard medical diagnostics**. Memorial Sloan Kettering describes bioresonance therapy as an **unproved** device approach used for diagnosis/treatment claims, and the FTC has acted against unsupported Zapper-related serious disease claims. Memorial Sloan ... +1

1. Zapper Device Datasheet

General Description

Item	Specification
Device category	Low-voltage electrical frequency generator
Common names	Zapper, Hulda Clark Zapper, parasite zapper, frequency zapper
Typical application field	Alternative wellness / frequency device
User interface	Usually switch, LED indicator, frequency selector, sometimes LCD
Contact method	Hand electrodes, wrist straps, foot pads, or alligator leads
Typical waveform	Positive-offset square wave
Typical frequency	30 kHz common; some units also include 1 kHz; adjustable models may cover wider ranges
Power source	Usually 9 V battery; some professional models use rechargeable or DC adapter
Output voltage	Common simple units: up to battery voltage, usually around 5–9 V depending on circuit and load
Output current	Usually low-current; exact value depends on model and body/contact impedance

Typical Electrical Specification

Parameter	Typical Value
Waveform	Square wave, often positive-offset
Frequency range	10 Hz–500 kHz claimed by some Clark-style definitions; common fixed models use 30 kHz and/or 1 kHz
Duty cycle	Commonly near 50%, model-dependent
Output coupling	Direct contact via electrodes/leads
Output terminals	Banana plug, alligator clip, copper hand electrodes, wrist straps
Battery	9 V alkaline battery
Enclosure	Plastic/ABS handheld enclosure

Parameter	Typical Value
Safety features	Battery indicator, insulated leads, low-voltage output; varies by manufacturer

A Clark-style Zapper vendor describes the device as producing a positive-offset frequency using a 9 V battery, with a dual-frequency model offering 30 kHz and 1 kHz, copper handles, banana-to-alligator leads, and a compact ABS enclosure. Clark Zapper

Mechanical Specification

Parameter	Typical Value
Size	Small handheld box
Example size	5.75 × 2.75 × 1 inch for one Clark-style model
Example weight	4.9 oz with battery for one Clark-style model
Accessories	Copper handles, cables, 9 V battery, pouch/case, instruction sheet

Typical Package Contents

Item	Included
Zapper main unit	Yes
Hand electrodes	Usually yes
Battery	Often yes
User manual	Usually yes
Carry pouch/case	Sometimes

Claimed Functions

Claimed Function	Notes
Frequency output	Real electrical function
Parasite/pathogen killing	Not clinically reliable; avoid medical claims
Detox support	Alternative/wellness claim

Claimed Function	Notes
General vitality support	Alternative/wellness claim

Safety / Compliance Notes

Area	Note
Implantable devices	Avoid use with pacemaker, ICD, implanted neurostimulator unless cleared by physician
Pregnancy	Avoid unless medically cleared
Broken skin	Avoid direct electrode contact
Serious disease claims	Do not use as treatment/diagnosis for cancer, infection, parasites, diabetes, etc.
Regulatory status	Many units are sold as wellness devices, not medical diagnostic equipment

2. MORA / MORA Nova Datasheet

General Description

Item	Specification
Device category	Bioresonance therapy / electrodermal measurement system
Representative model	MORA Nova / MORA Nova SRM
Manufacturer family	Med-Tronik / MORA system family
Application field	Bioresonance practice, electrodermal resistance measurement, alternative medicine
Use environment	Clinic/practice setting
Contact method	Electrodes, measuring stylus, foot electrode/console, accessory jacks
Display	Integrated touchscreen/display, depending on model
Measurement type	Skin electrical resistance measurement using very small voltage/current

Med-Tronik describes MORA Nova SRM as a device for measuring skin electrical resistance by applying a very weak electrical current via an electrode and counter-electrode.

Med-Tronik

Technical Specification — MORA Nova

Parameter	Specification
Dimensions	Width 470 mm × Height 147 mm × Depth 425 mm
Weight	Approx. 9 kg
Supply voltage	110–240 V AC, 50–60 Hz
Current input	700–300 mA
Buffer battery	Lithium battery BR2032
Classification	Class IIa, Type BF
Measurement range	< 820 Ω to 4.7 MΩ
Measurement accuracy	±1 scale unit
Measurement graduation	0–100
Max measurement current	< 25 µA DC
Max measurement voltage	< 1.25 V DC
Number of channels	2
Operating humidity	0–90%
Air pressure	60–106 kPa
Operating temperature	0–40 °C

These technical values come from the MORA Nova user manual technical-data section.

Med-Tronik

Interfaces / Connections

Interface	Purpose
Power socket	Mains power input
On/off switch	Main power
Measurement stylus jack	Electrodermal point measurement

Interface	Purpose
Foot console/electrode connection	Counter-electrode/contact measurement
Accessory jacks	Additional MORA accessories/modules
Touchscreen	Operation and measurement display

Measurement Principle

Item	Description
Measurement target	Skin resistance / electrodermal measurement point
Signal applied	Small measurement voltage between two electrodes
Measurement current	Very small current, limited by device
Output result	Numeric value, typically displayed on 0–100 scale
Interpretation	Practitioner-dependent bioresonance/electrodermal interpretation

The user manual states that during resistance measurement the device applies a small measurement voltage between two electrodes and measures the very small current, converting it into a displayed number.

Med-Tronik

Basic Accessories

Accessory	Purpose
Power cord	Mains operation
Foot electrode/console	Patient contact/counter-electrode
Measurement stylus	Point measurement
Touch pen	Screen operation
Accessory cables	Module/electrode connection
Trolley/case	Transport, depending on package

Claimed Functions

Claimed Function	Notes
Bioresonance analysis	Alternative-medicine interpretation

Claimed Function	Notes
Electrodermal resistance measurement	Actual electrical measurement
Therapy/program modules	Model and software dependent
Allergy/chronic condition support	Claimed by practitioners; not equivalent to validated clinical diagnosis

Safety / Maintenance

Area	Specification
Safety checks	Manual lists safety-related checks every 24 months
EMC	Manual includes IEC 60601 EMC guidance
User group	Manual warning says the system should be used by medical professionals
Cleaning/accessory handling	Follow manufacturer manual

The MORA manual includes EMC guidance and states the device/system must only be used by medical professionals.

Med-Tronik

3. Metatron / Metadron / Hunter 4025 NLS Datasheet

General Description

Item	Specification
Device category	NLS / non-linear analysis system
Common names	Metatron, Meta Hunter, Hunter 4025, Metatron 4025 Hunter
Typical application field	Alternative diagnostic-style body scanning / wellness analysis
Main components	Hardware interface box, headphones/bioinductor, software, laptop/PC
Contact method	Usually headphones/bioinductor placed on head
Software functions	Body atlas, organ/system selection, entropy analysis, etalon comparison, reports

Item	Specification
Claimed measurement type	"Wave characteristics," "spectral similarity," informational/energetic analysis

The Metatron Hunter user manual describes NLS systems as tracing body conditions through changes in wave characteristics of tissues/cells and includes software sections such as catalogue, etalons, entropy analysis, NLS analysis, meta-therapy, and re-printer.

ISHA Bioresonan... +1

Typical Hardware Specification

Parameter	Typical / Common Value
Main unit	USB-connected NLS interface box
Sensor/output device	Headphones / bioinductor
Computer	Windows laptop or desktop, model-dependent
Connection	USB, sometimes remote/software dongle
Power	USB-powered or external adapter, depending on unit
Display	Computer software interface
Data storage	Software database/patient card index
Report output	Software-generated report, usually printable/exportable

Typical Software Modules

Module	Purpose
Patient card index	Stores client records
Anatomical atlas	Selects organs/body systems
Catalogue	Database of reference patterns
Etalons	Reference items such as minerals, remedies, organs, pathologies, etc.
Entropy analysis	Claimed mathematical similarity/stage analysis
NLS analysis	Main scan/interpretation module
Meta-therapy	Claimed frequency/informational output

Module	Purpose
Vegeto test	Simulated compatibility/influence testing
Reprinter	Claimed transfer of information/frequency to carrier

The manual specifically shows menu sections for "Etalons," "Entropy analysis," "NLS Analysis," "Meta-Therapy," "Vegesto Test," and "Reprinter."

ISHA Bioresonan...

Claimed Functions

Claimed Function	Description
Whole-body scan	Claimed organ/system status visualization
Spectral comparison	Compares client scan to internal database/etalons
Entropy analysis	Claims to model pathological process stages
Remedy/etalon matching	Claims to compare minerals, remedies, microorganisms, tissues, etc.
Meta-therapy	Claims to transmit corrective/informational signal

The Metatron manual describes use of a bioinductor/headphones during "Vegeto-test" and discusses etalons and entropy analysis as part of its software workflow.

ISHA Bioresonan... +1

Typical Package Contents

Item	Included
NLS main box	Yes
Headphones / bioinductor	Yes
Software	Yes
USB cable	Yes
Laptop	Sometimes included by seller
Dongle/license key	Sometimes
Manual/training videos	Often
Carry case	Sometimes

Technical Gaps / Caution

Area	Note
Exact electrical output	Often not clearly specified in public brochures
Sensor principle	Usually described in vendor-specific non-mainstream terms
Diagnostic validation	Not equivalent to ultrasound, ECG, CT, MRI, blood testing, or FDA-cleared diagnostics
Supplier quality	Market includes many clones/copies; request serial number, manufacturer, software version, CE/FDA documentation if purchasing

4. Rayometer / Rayonex / Rayocomp PS 1000 Datasheet

General Description

Item	Specification
Device category	Bioresonance device according to Paul Schmidt
Representative model	Rayocomp PS 1000 polar 4.0
Related names	Rayometer, Rayonex, Rayocomp, Rayoscan, Rayotensor
Application field	Professional bioresonance/wellness practice
Use environment	Stationary practice/clinic use
Core concept	Frequency spectra / energetic analysis / harmonisation
Expansion	Modular software and accessories

Rayonex describes the Rayocomp PS 1000 polar 4.0 as a high-end bioresonance device for stationary practice use and says it can be expanded with modules such as RAH Human and Rayoscan Human. rayonex.co.uk

Product / Model Options

Model / Package	Description
Rayocomp PS 1000 polar 4.0 base unit	Main stationary device
Rayocomp PS 1000 polar 4.0 med.	Medical/practice version marketed as Class IIa in the UK listing

Model / Package	Description
Complete Premium package	Includes main device, software modules, chair/table, detectors, accessories, PS 10 portable device
Rayocomp PS 10	Portable/home-use companion device

A Rayonex UK product page lists the Rayocomp PS 1000 polar 4.0 Med as a Class IIa medical device with integrated RAH Human module and modular expansion capability.

rayonex.co.uk

Functional Specification

Function	Description
Frequency testing	Identifies claimed energetic frequency ranges
Harmonisation	Transmits frequency spectra via detectors
RAH Human	Human frequency/program database
Rayoscan Human	Analysis module using detector/ECG pad connection
Acupuncture oscillation therapy	Optional module/program set
PDF/printing	Some packages support PDF generation/printing
Language modules	German, English, French, Italian, Spanish, Dutch, Chinese, Japanese, Russian, Hebrew listed by Rayonex

Rayonex states that RAH Human includes over 1,800 frequency spectra for energetic testing/harmonisation and that Rayoscan uses connection via detector and four ECG pads to evaluate qualitative heartbeat changes. rayonex.co.uk

Hardware / System Specification

Parameter	Specification
Main device type	Stationary high-end bioresonance device
Power	External power supply/charger
Battery	Integrated/rechargeable battery used as short-term UPS
Connectivity	Wi-Fi USB adapter in some packages
Output/contact	Detectors, fabric detector, ECG pads, cups, Rayotensor

Parameter	Specification
Printing/export	PCL6 printer support and PDF generation described in Rayonex material
Manufacturing	Rayonex states devices are designed and manufactured in Lennestadt, Germany
CE	Rayonex states its devices carry CE marking

Rayonex lists an integrated battery as a UPS for short power interruptions, PDF generation/printing support, language modules, Made in Germany, and CE marking statements. [Rayonex Biome...](#)

Example Package Contents — Rayocomp PS 1000 polar 4.0 2025 Base

Item	Included
Rayocomp PS 1000 polar 4.0	Yes
Integrated polarisor	Yes
Manual / introduction book	Yes
Battery	Yes
Wi-Fi USB adapter	Yes
Protection cup / measuring cup	Yes
USB stick for PDF printing	Yes
Touch pen	Yes
Power supply/charger	Yes
Membrane packaging	Yes

The 2025 base-unit listing includes the main Rayocomp PS 1000 polar 4.0, integrated polarisor, manual, battery, Wi-Fi USB adapter, stainless cups, touch pen, power supply/charger, and related accessories. [rayonex.co.uk](#)

Example Package Contents — Complete Premium

Item	Included
Rayocomp PS 1000 polar 4.0	Yes
RAH Human	Yes

Item	Included
Rayoscan Human	Yes
RAH Global Test Human	Yes
RAH C-Module	Yes
Acupuncture oscillation therapy	Yes
Rayoscan cable and ECG holders	Yes
Fabric detector sheet	Yes
Fabric detector straps	Yes
Treatment chair/table	Yes
Rayocomp PS 10 basic	Yes
Green Cards	Yes
Rayotensor	Yes
Protection/measuring cups	Yes
Trolley/cases	Yes

Rayonex’s Complete Premium listing includes the Rayocomp PS 1000 polar 4.0, RAH/Rayoscan modules, detectors, Green Cards, Rayotensor, PS 10, treatment furniture, and accessories. [Rayonex Biome...](#)

Price Reference

Model	Listed Price
Rayocomp PS 1000 polar 4.0 2025 base unit	£21,990 on a UK listing

The UK listing shows £21,990 for the Rayocomp PS 1000 polar 4.0 2025 Edition base unit. [rayonex.co.uk](#)

5. NLS / Non-Linear Health System Datasheet

General Description

Item	Specification
Device category	Non-linear system / non-linear screening / NLS analyzer

Item	Specification
Common examples	Metatron Hunter, Meta Hunter, 3D/8D/18D NLS analyzers
Application field	Alternative wellness scan / energetic analysis
Main components	PC software, USB hardware box, sensor/headphones
Output	Body-system charts, organ images, indexes, risk-style reports
Claimed basis	Non-linear analysis, entropy, spectral database matching
Regulatory status	Usually not validated as mainstream diagnostic equipment

Typical System Architecture

Component	Description
Client software	Runs analysis workflow and database
Hardware interface	Connects PC to headset/sensor
Headphones/bioinductor	Used as sensor/output component
Database	Internal etalons/reference patterns
Report module	Exports or prints scan results
Optional therapy module	Claimed frequency/informational output

Typical Software Features

Feature	Description
Whole body scan	Selects body systems and organs
Organ visualization	Uses anatomical images/models
Entropy score	Claimed stage/severity indicator
Etalon comparison	Compares to remedies, diseases, tissues, microorganisms, etc.
Report generation	Creates client report
Multi-language support	Depends on seller/software

Feature	Description
Remote scan claims	Some sellers claim remote/off-site functions; treat cautiously

Typical Hardware Specification

Parameter	Typical Value
PC operating system	Commonly Windows, model/software-dependent
Connection	USB
Sensor	Headphone/bioinductor style
Power	USB or external DC adapter
Enclosure	Small desktop box
Accessories	Headset, USB cable, software/license, manual

Technical Limitations

Limitation	Detail
Not a standard diagnostic modality	Does not replace lab tests, ECG, X-ray, ultrasound, CT, MRI
Public specs often incomplete	Many sellers provide marketing claims, not engineering specs
High clone risk	Same names may be sold by many suppliers with different hardware/software
Accuracy claims	Often vendor-stated and not equivalent to independent clinical validation

6. Comparison Table

Feature	Zapper	MORA Nova	Metatron/NLS	Rayocomp/Rayonex	Generic System
Main type	Electrical pulse generator	Bioresonance + electrodermal measurement	Software-based NLS scanner	Bioresonance/harm onisation system	NLS sca categor

Feature	Zapper	MORA Nova	Metatron/NLS	Rayocomp/Rayonex	Generic System
Contact method	Hand/foot electrodes	Stylus + electrode	Headphones/bioinductor	Detectors, pads, cups, Rayotensor	Headpl
Main output	Electrical square wave	Resistance value + therapy modules	Software analysis report	Frequency/energetic test + harmonisation	Softwar report
Typical power	9 V battery	110–240 V AC	USB/PC	Power supply + battery UPS	USB/PC
Clinical diagnostic status	Not accepted	Alternative/electrodermal	Alternative/NLS	Alternative/bioresonance	Alterna
Engineering transparency	Simple	Better documented	Often poor	Moderate/vendor-specific	Often p
Best use classification	Wellness gadget	Professional alternative device	Wellness/NLS analysis	Professional bioresonance system	Wellnes softwar

7. What You Should Ask Supplier Before Buying

Use this checklist. It is very important because many NLS/bioresonance devices have unclear model names.

Required Document	Why It Matters
Exact model name and version	“Metatron” or “NLS” alone is not enough
Manufacturer name and address	Needed for traceability
Serial number format	Helps detect clone/copy products
User manual	Operation and safety
Technical datasheet	Real electrical/mechanical specs
Software version	NLS devices depend heavily on software
Certificate of conformity / CE declaration	Regulatory traceability
Medical device class claim	Must match region and intended use

Required Document	Why It Matters
Warranty terms	Many devices are expensive
Training material	These systems are operator-dependent
Output voltage/current limits	Important for electrical safety
Contraindications	Pacemaker, pregnancy, epilepsy, implanted devices, etc.
Service/calibration process	Especially for MORA/Rayonex-style systems

8. My Practical Recommendation

For **technical evaluation**, the safest ranking is:

Rank	Device family	Reason
1	MORA Nova	Has clearer technical data: dimensions, voltage, measurement range, current, voltage, channels
2	Rayocomp/Rayonex	More professional ecosystem and package information, but many specs are still vendor/module-dependent
3	Zapper	Electrically simple and easy to inspect, but health claims are weak/unsupported
4	Metatron/NLS	Software-heavy, many clone products, limited public engineering details
5	Generic non-linear health system	Too broad; must identify exact model first

For **medical use**, do not treat any of them as a replacement for approved diagnostic equipment. FDA has warned against certain unreviewed wellness/energy-style devices when safety and effectiveness were not established for diagnosing or treating medical conditions. [fda.gov](#)

 Sources

generate the picture for each device



could you give me the frequency list for disease which is used in zapper?

I can give you a **claimed / historical zapper frequency list**, but I'll keep it clearly separated from real medical facts. These frequencies are not clinically proven to diagnose or treat disease, so they should be treated as alternative-device reference data only.

Thought for 10s >

Yes — but important distinction:

Classic Hulda Clark “zapper” usually uses a **fixed frequency**, commonly **30 kHz** square wave, sometimes with **1 kHz** as a second mode. Disease-by-disease frequency lists are more commonly from **Rife / CAFL / multi-frequency zapper** systems, not the simple original zapper. Some Clark-style sources also describe scanning ranges such as **290,000–470,000 Hz**, but these are alternative-device claims, not validated medical treatment protocols.

Hulda Clark Para... +1

Also, these lists are **claimed/experimental** frequency lists. They are **not clinically proven** to diagnose, cure, or treat disease. The FTC has taken action against unsupported zapper-related disease claims, so please use this only as technical/reference information, not medical advice.

Fede... +1

Common Zapper / Rife-style Frequency Reference List

General zapper frequencies

Use category	Claimed frequency
Classic Hulda Clark zapper	30 kHz
Some dual-frequency zapper devices	30 kHz / 1 kHz
Clark-style multi-frequency sweep	290 kHz–470 kHz
Some broader frequency generators	Hz to kHz/MHz range, model-dependent

Disease / condition frequency list — claimed reference only

The following values are examples from public zapper/Rife-style frequency lists. Many lists give frequencies in **kHz**, so multiply by **1000** to convert to Hz. For example, **7.50 kHz = 7,500 Hz**. One public “Zapper Frequency List” explicitly says its list is in **kHz**. Detoks

Condition / disease name in list	Claimed frequencies
Addison disease	0.19, 0.75, 0.90, 7.50, 27.50, 222.70, 425.71, 563.19, 642.91, 978.05 kHz
Adenoma	0.04, 5.81, 22.50, 52.50, 92.50, 224.37, 527.00, 667.00, 721.00, 987.23 kHz
Abdominal cramps / stomach	72, 95, 190, 304 Hz in one Rife-style list
Lactobacillus acidophilus	865, 17375.68 Hz in one Rife-style list
Inflammation, general	Often listed differently by different sources; no single standard
Parasites, general	Often zapper sellers use 30 kHz or frequency sweeps rather than one exact frequency
Bacteria, general	Often listed under organism-specific names, not one universal frequency
Virus, general	Usually listed by specific virus name in Rife/CAFL lists
Candida / yeast	Frequency varies widely by list
Lyme / Borrelia	Frequency varies widely by list
Staphylococcus	Frequency varies widely by list
Streptococcus	Frequency varies widely by list
E. coli	Frequency varies widely by list

Condition / disease name in list	Claimed frequencies
Herpes	Frequency varies widely by list
HPV	Frequency varies widely by list
Influenza	Frequency varies widely by list
Sinusitis	Frequency varies widely by list
Arthritis	Frequency varies widely by list
Acne	Frequency varies widely by list
Migraine / headache	Frequency varies widely by list

Why I cannot give one “official” disease frequency table

There is **no medically accepted official disease-frequency list** for zapper use. Public lists come from different alternative sources such as Clark-style lists, Rife lists, CAFL lists, Spooky2-style databases, and vendor-specific programs. The CAFL page itself frames the list as experimental, and Rife-frequency sites usually describe these as research or alternative frequency references rather than approved medical protocols.

Electroherbalism +1

Best practical format

For a real zapper/frequency device datasheet, I recommend structuring the frequency table like this:

Field	Example
Program name	Parasite general
Frequency unit	Hz or kHz
Frequency list	30,000 Hz or sweep 290,000–470,000 Hz
Waveform	Square wave
Offset	Positive offset, if Clark-style
Voltage	Model-dependent
Current limit	Model-dependent
Source	Clark / Rife / CAFL / vendor
Medical status	Not clinically validated

The most honest summary: simple zappers usually use 30 kHz, while disease-specific frequency lists belong more to Rife/CAFL-style frequency generators.

      Sources

then generate the picture for the signal which is generated in the zapper.

Stopped thinking



   