







Modules

Since its first development, Mora NOVA has been continuously improved. The fact that the computer program is built in the device allows the modules to be easily installed and used on the device. Module options can be purchased additionally in line with users' needs and can be easily installed on user devices via the internet. Some modules are offered to users as standard in the delivery of the Turkish Standard Device. These are;

- Diagnostic Module
- Optima Module
- Point Therapy Module
- Universal Module
- Tensor Module

Additionally, modules that can be purchased according to their intended use and are almost not found in any other bioresonance device;

- Cornelissen Module
- 4-Pol Module
- Psychobiophonie Module
- Dental Module
- 5 Elements Music Module
- Element Group Module

Technical Specifications

- Dual Channel Technology, 2 Interfaces (A Mode and A-bar Mode),
- Significant improvement in the Scott-Morley test for 2 Channels,
- Frequency range: 0.1 HZ to 1MHz, Filter Adjustment Range: 1 Hz to 400 kHz,
- Amplification between 0.1 and 1 Million for each Channel and Mode,
- Modular Design, Computer is inside the Device,
- Programmable feature up to 16 individual steps,
- Ability to automatically measure from 4 to 7 segments,
- Automatic detection of acute and chronic conditions by the device,
- Patient Registration System,
- Integration of Mora Mouse into the device,
- Freely adjustable therapy from 1 to 65000,
- Automatic display of Measurement Points graphically on the screen,
- Laser Electrodes,
- Square Wave Generator 1 Hz to 500 kHz, Sine Wave Generator 1 Hz to 250 kHz.,
- 15 inch LED illuminated touch screen,
- Printable feature.

Mora Nova

The latest product of bioresonance technology; The Perfect Combination of Design and Technology

Mora Therapy has taken the successful point it achieved with Mora Super Plus Devices significantly further with MORA NOVA devices. As you know, every Technology requires innovation and change. In bioresonance Technologies, a technological development as important as a revolution has been achieved with Mora Nova. The usage has been greatly facilitated with its modern touch screen interface, which is unmatched in any counterpart. In addition to its built-in computerized touch screen, the technical specifications of the bioresonance device used inside have been greatly improved, leading to a much stronger treatment effectiveness. It contains many different and comprehensive test sets, modules, therapy & diagnostic kits, and also brings many innovations not only to bioresonance but also to holistic medical technologies.

Striking Advantages of Mora NOVA;

Mora NOVA is the latest and innovative version of the original bioresonance therapies developed by Dr. Franz Morell and Erich Rasche, utilizing the most advanced and up-to-date technology both in device and software. Mora NOVA has a respected position in the world with its advanced technology applications, ease of use and high efficacy in Bioerzonance devices.

- Mora Nova is much different and more effective than Mora SP and other biorasonance devices in terms of frequency scanning range. This superiority provides higher therapy success rates in all applications with the new Mora Nova.
- It works through a dual-channel software program that operates on 16 separate tabs, which allows to apply both color therapies, test kits and programs at the same time.
- It allows automatic operation in a wide range of different indications with the "Mora Optima" module, which offers new and functional user-specific therapies.

- All patient information can be easily accessed through the data storage feature in Mora NOVA when needed.
- The new program interface provides the user with the opportunity to perform allergy therapies and tests much more easily. All other new tests such as food intolerance, vega, amalgam, candida and heavy metal etc. can be applied much more easily.
- In the Turkish Standard set; Mora Color, Bach Flowers, System Analysis Test Kit, Optima, Point Therapy, Universal Therapy modules, Diagnostic Module are offered to users in addition to classical bioresonance treatments.
- With the Basic Therapy Module, also known as Therapy Preparation, there is a program interface that offers special therapy opportunities according to the age range and dominant problem of the client.
- During therapy, the patient's acute and chronic condition can be automatically selected by New Mora Nova.
- Thanks to its BlueMag Magnetic Electrode, it provides high efficiency and ease of use in local applications with therapies.
- Without the need for additional equipment, it can easily use Color Therapy, Bach Flowers Therapy, and all other test kits with its built-in computer.
- Thanks to its innovative and modern design, it has a 15-inch foldable, touch-sensitive, high-resolution LED screen that provides easy access to all menus and programs.
- The body part undergoing Diagnosis and Therapy can be seen on the screen as graphics and images, and the points at which the application can be performed are automatically displayed via the guidance menu.
- With the graphical menu navigation, all applications can be easily accessed from the main screen
- Due to its one-piece compact design, it takes up less space and the device is very easy to carry.
- It has a practical use feature with very little use of cables - The cups are fixedly placed in their special slots in the front part of Mora Nova, thus providing a more reliable usage feature.
- Foot Electrode can be mounted on the device during transportation or waiting.
- It can be directly connected to the Internet, so remote assistance, updates and upgrades can be easily performed.

- With the company's specially developed Satellite program for the device, it can be connected to a computer or laptop, and all information in Mora Nova can be printed via this program, and therapy before and after comparisons can be made with graphic programs.

Test Sets

With Mora NOVA, some test sets have been studied extensively, new items have been added to some of the test sets, some tests have been combined to make them much easier to use, and new test sets that were not previously included in Mora applications have been added.

In general, we can list the Test Sets as follows;

- Nova Elements Group Test Set: 5 Elements, Allergy, Chemical Elements, Geopathy, Mycosis, Nosodes, Sussler Salt, Horvi, Meridian Test Sets.
- Homeo Group Test Set: Heel, Wala, Vitortho, Orthica, Inus Test Sets.
- Ortho Group: Orthica, Inus, Vithorto, Sanum Test Sets.
- Cornelissen: System Diagnostics, Mitochondropathy, Cancer, Alzheimer, Dementia Test Sets.
- AllerHomeo Group: Allergy, Wala Test Sets
- Allergy Test Set; There have been comprehensive changes in the Allergy Test Set in the New Mora NOVA. Although the General Structure seems the same, many new items have been added; and the Dermatology and Pesticides Test Sets, which were previously separate, have been combined
- Bach Flower Test set
- Color Therapy Test Set
- Dr. Cornellissen System Diagnostic Test Set
- Dr. Cornellissen Mitochondropathy / Cancer Test Set
- Dr. Cornellissen Alzheimers + Dementia Test Set
- Heel Test Set; This Test kit contains all Heel Homeopathic medicines and Biological Catalysts.
- Geopathy Test Set
- Homeopathic Potencies and Dilutions - Korsakov Test Set; A whole has been created by adding Hautes Dynamisations and Homeopathic Substances Test Set to this Test Set.
- 5-Element Test Set
- Meridian Test Set

- Schüssler Salts Test Set
- Nosodes Test Set
- Trace Elements Test Set; This Test Set also contains both trace elements and chemical elements.
- Mycosis Test Set
- Horvi Test Set
- Wala Test Set
- Orthica Test Set
- Inus Test Set
- Vitortho Test Set
- Sanum Test Set
- Traditional Chinese Medicine Test Set
- Pascoee Test Set

Share



EDT Sağlık Hizmetleri Tic. Ltd. is the sole authorized distributor of MORA Technology in Turkey and Northern Cyprus.



Contact

Fahrettin Altay Mah. Mithat
Paşa Cad. No:1163 Halil Bey
Apt Kat:2 Daire:4 35140
Karabağlar / İzmir

Working hours

Monday-Friday:
09:00-18:00
Saturday Sunday:
Closed

✉ info@mora.com.tr

☎ +90 216 405 14 52 - +90 216
405 14 53

☎ +90 533 814 50 11



Mora Therapy in the media

2023 © EDT Sağlık Hizmetleri Tic. Ltd. Şti. All Rights Reserved. Web Design: Quattro Project
Personal Data Processing | KVKK Information Text | KVKK Application Form